## FR – Instructions d'installation

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**EN – Installation instructions**

---

You will find below the link and the instructions to download and install the Global D libraries. For more information, visit our website <https://www.globald.com/en/implantology/implantology-digital-solutions/cad/>

V3 libraries are validated by exocad

**Library download link from the exocad website::**

- Link : <https://exocad.com/integration/dentalcad-library-integration#/all/tag-g>
- Select packages (+ Add to download list)
- GlobalD MultiUnit exocad libraries (10 libraries)
- GlobalD IL exocad libraries (2 libraries)
- GlobalD Lab Analog exocad libraries (1 library)

It's possible to maintain the old versions V2.2e & V2.4.4 (Global D V2) but it is preferable to delete them.

or

**Library download link from Global D**

Link : <https://monextra.net/wl/?id=8BMhxEF14rggBqCbuKBk5jixGjZOMxDo>

Mot de passe : Gld2018-Exocad

Libraries & documentation

**Find the installation folders location:**

The 2 locations to find are:

- ...\\DentalCADApp\\library\\implant\\ for single or plural abutment libraries,
- ...\\DentalCADApp\\library\\modelcreator\\implant\\ for the model library.
- The link to the DentalCADApp folder location can be identified by right-clicking on the Exocad executable icon (cad), then go to Properties and noting the path in the "Target" box, such as "C:\\Exocad\\exocad-DentalCAD....\\DentalCADApp\\...". This location may be different with some integrators, for example:
- Amann Girrbach : ...\\Ceramill software\\CAD\\library\\
- Zfx : ...\\Zfx Dental Design\\ZfxCAD\\library\\

**Exocad libraries installation :**

*do not change folders names*

- Unzip all files
  - Go to ...\\DentalCADApp\\library\\implant\\
-

- Delete older Global D libraries if necessary (optional),
- Drop folders

For CAD	For Exoplan
Global-D_EVL_conid	Global-D_twinKon_conid_multiunit
Global-D_In-Kone_conid	Global-D_EVL4_conid_multiunit
Global-D_MUA_conid_multiunit_support	Global-D_EVL5_conid_multiunit
Global-D_MULTI_conid_multiunit_support	Global-D_EVL33_conid_multiunit
Global-D_PAOF_conid_multiunit_support	Global-D_In-Kone_ST_conid_multiunit
Global-D_PC43_conid_multiunit_support	Global-D_In-Kone_WD_conid_multiunit

- Go to ...\\DentalCADApp\\library\\modelcreator\\implant
- Delete older Global D libraries, if necessary (optional),
- Déposer le dossier "Global-D\_ModelCreator"

All folders must contain files (.fda / .stl / .xml / .txt) and no folder. If there is one or more folders.

#### Selection Menus

- "Global-D - xxxx" New libraries (V3) to use
- "Global D V2" if maintain, old complete V2.2e library
- "Global D WD" if maintain, add-on with only the prosthetic connection on In-Kone WD implant

#### The CAD-CAM connections provided:

- **In-Kone ST** on implant and indexed Ti Base In-Kone Ø3.8 Gingival heights 0.6, 2, 3 and 4 mm or Ø5.5 heights 1.5, 2.2, 3 et 4 mm,
- **In-Kone WD** on implant and indexed Ti Base In-Kone for healing screws Ø6.5, Gingival heights 1.5, 2.2, 3 et 4 mm,
- **EVL** on implant and indexed Ti Base EVL (except EVL C),
- **EVL NI** on implant and Ti Base non-indexed EVL (except Exocad – 100% panthera)
- **PC 4.3** on new conical abutments 0 ° and rights associated sockets, non-indexed (In-Kone, TwinKon)
- **PAOF** on angulated conical abutments and associated sockets, non-indexed In-Kone (In-Kone, compatible with Hexacolor Ø4.7 angulated conical abutment)
- **MULTI** on MULTI abutment and bases, implant EVLC, non-indexed (EVL, CLINIC, compatible with In-kone PC4.0 and Hexacolor conical abutment)
- **MUA** on pillars and sockets MULTI angulated, non-indexed (EVL, CLINIC,).
- **In-Kone+PC4.3 or PC4.3 implant Studio** reserved in synergy with a guided surgery planning,
- **EVL + MULTI** reserved in synergy with a guided surgery planning.

#### Examples of selection:

- For designing an abutment on titanium base, please select ...> **BaseTi**
- For designing an abutment on titanium base 2 mm shortened on lab (top shaft = 3 mm instead of 5 mm (standard)), please select ... > **BaseTi Ø... (Hc 3mm)**
- For designing a customized abutment from premilled on S3DEL assembly on lab (In-Kone example), please select **Global-D - In-Kone > In-Kone ST > Premilled S3DEL**
- For designing a customized Bridge or bar « Compliance » without Ti Base (PC 4.3 example), please select **Global-D - PC4.3 > Compliance**

For more information see file "DigitalSol\_CAD-CAM\_QuickStart..."

This library includes several heights of In-Kone bases and offers the possibility to the lab to shorten them to 3 mm (i.e. at the level of the middle groove). **Warning 3D direct connections implants or abutments are false. Only the Panthera Dental Company is authorized to machine these connections. Machining Laboratory is possible via Tibase, premilled or compliance system.**

We count on you to let us know your suggestions for correction or improvement.