FR – Instructions d'installation

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EN – Installation instructions

3Shape folder download link containing the libraries:

- Global_D_3Shape_library_Complete_V3.0.9_Pack.zip
[Public link to the library](#)
(Contains dme files, multi-language documentations and installation instructions)

Unzip the following files:

- Global_D_3Shape_library_Complete_V3.0.9_2020.dme (3Shape Dental System 2020)
- Global_D_3Shape_library_Complete_V3.0.9_2021.dme (3Shape Dental System 2021)
- Global_D_3Shape_library_Complete_V3.0.9_2022.dme (3Shape Dental System 2022 or more)

First install :

- Run "Dental System Control Panel",
- Go to "Tools" and select "Import/Export"
- Click to "import materials"
- Select the file "Global_D_3Shape_library_Complete_V3.0.9_2022.dme"
- Return to the home page
- Go to "Site Setting" and select "Manufacturing processes"
- Select (left table) "Hard machining" and/or "Milling R0.4" and/or "Milling R0.5" if they exist, or the usual process that you select for metal machining
- Associate the manufacturer (you in right table) with the Global D Materials (center table) by checking "Global D Ti Premilled", "Global D CoCr Compliance" and "Global D Ti Compliance"
- Exit the System Control Panel and "save changes before closing"

Update of 3Shape library:

- Run Dental "System Control Panel"
 - Go to "Abutments" and select "implant system"
 - Check Global D library versions (in the list starting from the end)
-

-
- Delete all Global D implant system, version greater to V3.0.9
 - Return to the home page
 - Go to "Basic Elements" and select "Materials",
 - Select and delete in the left table all the materials starting with Global D....
 - Save (disk in header menu),
 - Return to the home page
 - Go to "Tools" and select "Import/Export"
 - Click to "import materials"
 - Select the file "Global_D_3Shape_library_Complete_V3.0.9_2022.dme"
 - Return to the home page
 - Go to "Site Setting" and select "Manufacturing processes"
 - Select (left table) "Hard machining" and/or "Milling R0.4" and/or "Milling R0.5" if they exist, or the usual process that you select for metal machining
 - Associate the manufacturer (you in right table) with the Global D Materials (center table) by checking "Global D Ti Premilled", "Global D CoCr Compliance" and "Global D Ti Compliance"
 - Exit the System Control Panel and "save changes before closing"

Selection menus:

- "Global D" to be deleted,
- "Global D V2", to be deleted,
- "Global D WD" to be deleted,
- "Global D V3" contains V3.0.9 library to use

Warning

The "GLOBAL D In-Kone + PC4.3 " implant kit is only reserved for Implant Studio planning.

The CAD-CAM connections provided:

- In-Kone ST on implant
- In-Kone WD on implant,
- twinKon on implant,
- EVL on implant (without aesthetic ring and except EVL C),
- PC 4.3 on new conical abutments 0 °(In-Kone, TwinKon)
- PAOF on angulated conical abutments (In-Kone, compatible with Hexacolor Ø4.7 angulated conical abutment)
- MULTI on MULTI abutment and implant EVLC, (EVL, CLINIC, compatible with In-kone PC4.0 and Hexacolor conical abutment)
- MUA on angulated MULTI abutments (EVL, CLINIC,).
- In-Kone + PC4.3 reserved in synergy with guided surgery planning.

Examples of selection:

- For designing an abutment on titanium base, please select ...> ...**BaseTi**
- For designing an abutment on titanium base 2 mm shortened on lab (top shaft = 3 mm instead of 5 mm (standard)), please select ... > ...**BaseTi (Hc=3mm)**
- For designing a customized abutment from premilled on S3DEL assembly on lab (In-Kone example), please select **Global D V3 > Global D In-Kone > Premilled S3DEL**,
- For designing a customized Bridge or bar « Compliance » without Ti Base (PC 4.3 example), please select **Global D V3 > Global D PC4.3 > PC4.3 Compliance**

For more information see file "DigitalSol_CAD-CAM_QuickStart_..."

This library includes several heights of In-Kone bases and offers the possibility to the lab to shorten them to 3 mm (i.e. at the level of the middle groove). **Warning 3D direct connections implants or abutments are false. Only the Panthera Dental Company is authorized to machine these connections. Machining Laboratory is possible via Tibase, premilled or compliance system.**

We count on you to let us know your suggestions for correction or improvement.